



Email : sales@laboratoryinstrumentindia.com

Product Name :
Calibration Trainer

Product Code :
LBNY-0005-12700095



Description :

2 DC voltage sources to supply control loop components with auxiliary power
2 precision manometers
2 pressure regulators with manometers
Height-adjustable tank and scale to set low pressures
2 controllers to send and measure current signals
Investigation of transmission behaviour and calibration of control loop components and measuring instruments
Sending and measuring pneumatic and electrical signals
Control loop components and meters available as accessories.

Technical Specification :

2 pressure regulators
0...6bar
0...1,6bar
Pressure range, height-adjustable tank
0...1000mm head
1 output each: 4...20mA
1 input each: 4...20mA
2 controllers
2 DC voltage sources
each 24VDC
2 precision manometers
0...6bar, D=160mm

0...1,6bar, D=160mm.



Laboratory Instrument India