



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Maintenance Of Valves And Fittings And Actuators

**Product Code :**  
LBNY-0005-1270097



#### **Description :**

Pressure reducing valve  
Segmented ball valve with swivel drive  
2 compressed air ranges, adjustable by fine pressure regulator  
Electric signal transmitter for positioner in the form of an adjustable current source  
Maintenance work on industrial valves and fittings  
Pneumatic control valve with electro-pneumatic positioner  
Butterfly valve with swivel drive  
Instrumentation: analogue pressure meter, digital ammeter and voltmeter.

#### **Technical Specification :**

Single-action with spring return  
Pneumatic swivel drive  
Measuring ranges  
Pressure (bourdon tube manometer)  
0...2,5bar (D=250mm)  
0...6,0bar (D=60mm, fine pressure regulator)  
0...1,0bar (D=160mm)  
0...1,6bar (D=60mm, fine pressure regulator)  
Differential pressure: 0...10kPa  
Voltage (digital display): 0...20VDC  
Current (digital display): 0...20mA.



**Laboratory Instrument India**