



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Training System Flow Control Experiment

**Product Code :**  
LBNY-0005-127000103



**Description :**

Storage tank  
Capacity: approx. 3000mL  
Pump  
Max. head: 6m  
Max. flow rate: 8L/min  
Power consumption: 18W  
Proportional valve: Kvs: 0,7m<sup>3</sup>/h  
Flow sensor: 0,5...3L/min  
Rotameter: 20...250L/h  
Software controller configurable as P, PI, PID and switching controller.

**Technical Specification :**

Software  
Simulation function  
Disturbance variable input  
Process schematic with controller type selection  
(manual, continuous controller, two- or three-point controller, programmer)  
Time functions.

---

# Laboratory Instrument India

