



Email : sales@laboratoryinstrumentindia.com

Product Name :
Training System Flow Control Experiment

Product Code :
LBNY-0005-127000103



Description :

Storage tank
Capacity: approx. 3000mL
Pump
Max. head: 6m
Max. flow rate: 8L/min
Power consumption: 18W
Proportional valve: Kvs: 0,7m³/h
Flow sensor: 0,5...3L/min
Rotameter: 20...250L/h
Software controller configurable as P, PI, PID and switching controller.

Technical Specification :

Software
Simulation function
Disturbance variable input
Process schematic with controller type selection
(manual, continuous controller, two- or three-point controller, programmer)
Time functions.

Laboratory Instrument India

