




Email : sales@laboratoryinstrumentindia.com

Product Name :
Training System Position Control

Product Code :
LBNY-0005-127000105



Description :

Experimental unit for control engineering experiments
Rotary encoder as displacement sensor
Transparent protective cover
Linear position control of carriage with linear drive and gear motor
2 microswitches to shut down at end positions
Process schematic on front panel
Software-based controlled system simulation.

Technical Specification :

DC motor
Speed: 85min⁻¹
Torque: 200Nmm
Max. traverse rate: 45mm/s
Scale: 0...300mm
Transmission ratio: $i=50$
Travel: max. 300mm.

Laboratory Instrument India

