



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Smart Flow Transmitter Learning System

**Product Code :**  
LBNY-0005-1280006



**Description :**

The Smart Flow Transmitter Learning System includes an industrial quality transmitter and student learning materials for theory and lab. The transmitter mounts to the workstation console and connects to 3 optional flow transducers: pitot tube, venturi, and orifice plate. The Smart Flow Transmitter Learning System adds to Process Control Learning System to provide a 4-20 mA feedback signal that is proportional to the flow rate measured by transducers that sense differential pressure.

**Features:**

- Industrial-quality component
- Student Reference Guide
- Senses the differential pressure within the system.

**Technical Specification :**

Smart Flow Transmitter Learning System



**Laboratory Instrument India**

