



Email : sales@laboratoryinstrumentindia.com

Product Name :
Acidity (Ph) Control Apparatus

Product Code :
LBNY-0005-12800010



www.laboratoryinstrumentindia.com

Description :

pH Control Apparatus is designed to demonstrate the control of the acidity of a water supply by chemical dosing. Control of pH is an important aspect of water treatment in water cooling systems, hot water heating systems and boiler feed water. The pH Control Apparatus is designed to demonstrate the control of the acidity of a water supply by chemical dosing. Two methods are available to control the addition of the dosing chemical, either slug (on-off) dosing or continuous controlled injection. Both methods can be demonstrated with this apparatus. Optional equipment is available to convert the apparatus into a Ratio Control Apparatus. Two methods are available to control the addition of the dosing chemical, either slug (on-off) dosing or continuous controlled injection.

Technical Specification :

Acidity (Ph) Control Apparatus



Laboratory Instrument India