



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Optical Bench 3 Rails

**Product Code :**  
LBNY-0004-780031



**Description :**

Optical Bench 3 Rails

**Technical Specification :**

Optical bench 3 rails

Optical bench made of special aluminium profile, with excellent functionality, stability and low in weight. One long side with a metric mm scale, the other with a three-coloured cm block scale.

For computer aided experiments.

Length: 21 m

Metric scale in mm.



**Laboratory Instrument India**