



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Electronic Control Console

**Product Code :**  
LBNY-005-1280021



### Description :

Using this console also enables the sensor data to be logged and displayed on a computer (using the interface), but with the control functions being implemented in the hardware, rather than the computer software. The enables all of the single loop configurations to be implemented without using a computer.

### Technical Specification :

#### Facilities include:

- Use of 4-20mA interfaces.
- Use of serial communications (RS232) in a control system
- Simple interfacing to other PLC and/or PID controllers
- Use of a commercial PID controller
- Simple SCADA demonstration (requires PC).



**Laboratory Instrument India**