



Email : sales@laboratoryinstrumentindia.com

Product Name :
Statics Experiments Kit

Product Code :
LBNY-0005-1290003



Description :

All experiment elements are made in special anodized aluminium.

Centres of gravity of various shaped plates:

3 Cords and ring.

5 Cords and ring.

Beam balance.

Drawing panel.

Beam.

Pivot screw.

Adjustable hooks.

2 forces equality divided.

The simple pendulum.

Set of weights:

10 weights of 0.1 N.

2 weights of 0.5 N.

10 weights of 0.01 N.

4 weights of 0.05 N.

1 weight of 5 N.

4 weights of 1 N.

2 weights of 2 N.

Weight hooks.

Light weight hooks.

Screws.

Pulleys.

Large ext. spring.
Small ext. spring.
Knurled nuts.
Spare rope.
Dynamometer.

Technical Specification :

Statics Experiments Kit



Laboratory Instrument India