



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dynamics Experiments kit

Product Code :
LBNY-0005-1290004



Description :

This kit is useful for developing various models of physics. Advanced instruments of this kit enables one to understand flow characteristics of air and water, generated energy caused by their flowing and resting to name a few. Dynamics Experiments Kit is used for measuring various properties and effects of accelerating forces on inelastic and elastic collision, kinetic energy, masses and so on. We are an eminent supplier and manufacturer of Dynamics Experiments Kit. This kit comprises of instruments like force board, aluminum profiled trolley, cotton reel, ticker tape timer, spring balance and so on. This product has been designed to suit physics teaching needs of teachers and students as well.

Features:

Precise design of instruments
Sturdy construction of kit
Easy to handle.

Technical Specification :

Dynamics Experiments kit



Laboratory Instrument India

