




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Transmissions Experiments Kit

**Product Code :**  
LBNY-0005-1290005



**Description :**

This kit consists of suitable accessories for determining emission rate, color and ph level, software program manual and light source. Its application can be noticed in various teaching laboratories, earth science center and chemistry labs. Transmissions Experiments Kit is used for suitable measurement of different attributes of transmission lines and their input impedance. This multipurpose training system is used for providing thorough knowledge on phase delay, standing waves, attenuation and frequency response. This kit is suitable for explaining theories associated with color measurement based transmission, atomic emission and so on. We are a prominent manufacturer and exporter of Transmissions Experiments Kit.

**Features:**

Easy to understand mechanism  
Simple to use  
Meticulously designed.

**Technical Specification :**

Transmissions Experiments Kit



**Laboratory Instrument India**

