



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Load Elevation Mechanisms Experiments Kit

**Product Code :**  
LBNY-0005-1290007



**Description :**

All experiment elements are made in special anodized aluminium.

Wheel and axle set.

Weston differential chain block.

Three pulley block.

Support screw.

Two pulley block.

Screw jack.

Adjustable hooks.

Set of weights:

10 weights of 0.01 N

2 weights of 0.5 N.

4 weights of 0.05 N.

10 weights of 0.1 N.

4 weights of 1 N.

1 weight of 5 N.

2 weights of 2 N.

Weight hooks.

Light weight hook.

Single pulley block.

Knurled nuts.

Pulleys.

Spare rope.

---

Screws.  
Adjustable pulley.  
Dynamometer.

**Technical Specification :**

Load Elevation Mechanisms Experiments Kit



**Laboratory Instrument India**