



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Intermittent Motion Set

**Product Code :**  
LBNY-0005-1300007



**Description :**

The Intermittent Motion Set set is supplied in plastics box. A Geneva stop and a ratchet mechanism are used to demonstrate how rotary motion is converted into a start stop motion. Intermittent Motion Set is part of mechanical Engineering Test Series which provides sets for various experiments on Basic Panel (separately supplied) to study or demonstrate specific engineering principles.

**Technical Specification :**

**Experiments:**

Relation between rotary motion and start stop angular motion.



**Laboratory Instrument India**