




Email : sales@laboratoryinstrumentindia.com

Product Name :
Rotary Translation Motion Set

Product Code :
LBNY-0005-1300008



Description :

The set is supplied in plastics box. Slider crank, slotted link, and quick return mechanism are used for studying relation of rotary motion and translation. Rotary Translation Motion Set is part of mechanical Engineering Test Series which provides sets for various experiments on Basic Panel (separately supplied) to study or demonstrate specific engineering principles.

Technical Specification :

Experiments:

Relation of rotary motion and translation in slider crank, slotted link and quick return mechanism.



Laboratory Instrument India