



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bar Linkages Set

**Product Code :**  
LBNY-0005-1300009



**Description :**

The set is supplied in plastics box. Over 20 perforated bars of different lengths and joints are used for assembly to form different linkages. The bar can accommodate a pen to trace the movement of the linkages on paper. Bar Linkages Set is part of mechanical Engineering Test Series which provides sets for various experiments on Basic Panel (separately supplied) to study or demonstrate specific engineering principles.

**Technical Specification :**

**Experiments:**

Pantograph.  
Ackmann steering.  
Four Bar Linkages  
Straight line linkages.



**Laboratory Instrument India**