



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tensile Test Set

Product Code :
LBNY-0005-13000010



Description :

The set is supplied in plastics box. Tensile specimen is attached to the tensile test frame where load is a applied by a crew. This is part of Mechanical Engineering Test Series which provides sets for various experiments on Basic Panel (separately supplied) to study or demonstrate specific engineering principles. Load is measured by a dial indicator on the spring and specimen displacement is measured by a vernier.

Technical Specification :

Experiments:

Tensile test to failure to determine elastic limit, yield and ultimate strength of specimens.



Laboratory Instrument India