




Email : sales@laboratoryinstrumentindia.com

Product Name :
Crank And Toggle Set

Product Code :
LBNY-0005-1300013



Description :

The set is supplied in plastics box. Two parts beams with scale for bending moment and two parts beam with scale for shear force test with beam supports are used. Bending moment and shear force are calculated from spring balance readings, location of loads and weights. Crank And Toggle Set is part of Mechanical Engineering Test Series which provides sets for various experiments on Basic Panel (separately supplied) to study or demonstrate specific engineering principles. Loading is by load hangers and weights.

Technical Specification :

Experiments:

Shear force of a beam under loads.
Bending moment of a beam under loads.



Laboratory Instrument India