




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Potential And Kinetic Energy Set

**Product Code :**  
LBNY-0005-13000014



**Description :**

A pendulum, a flywheel with string, a weight hanger and weights are used to demonstrate how energy is transferred from one form to another. The set is supplied in a plastic box. Potential And Kinetic Energy Set is part of Mechanical Engineering Test Series which provides sets for various experiments on Basic Panel (separately supplied) to study or demonstrate specific engineering principles.

**Technical Specification :**

**Experiments:**

Static potential energy in a spring.  
Kinetic energy in a flywheel.  
Kinetic and potential energy in a pendulum.



**Laboratory Instrument India**