



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Optical Bench 3 Accessory 03

**Product Code :**  
LBNY-0004-780034



**Description :**

Optical Bench 3 Accessory 03

**Technical Specification :**

Optical bench 3 accessory 03:

It is made of smooth (fine polished) wooden base board with a wooden scale fixed on it. All accessories are mounted on a rod and can be accommodated in holders.

1 Lens (about +100 mm)

1 Lens (about +200 mm) (with holders if necessary)

1 Viewing screen.



**Laboratory Instrument India**