




Email : sales@laboratoryinstrumentindia.com

Product Name :
Spring Test Set

Product Code :
LBNY-0005-13000024



Description :

The set is supplied in plastics box. Coil spring(s) are attached to the test frame and load is applied to the spring(s). Spring extension is measured by a vernier. This is part of Mechanical Engineering Test Series which provides sets of various experiments on Basic Panel (separately supplied) to study or demonstrate specific engineering principles.

Technical Specification :

Experiments of Spring Test Set
Parallel and series spring tests.
Hooks law for spring tension.
Hooks law for spring compression.



Laboratory Instrument India