



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Optical Bench 3 Accessory 04

**Product Code :**  
LBNY-0004-780035



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Optical Bench 3 Accessory 04

**Technical Specification :**

Optical bench 3 accessory 04:  
Aperture wheel ? 6 slits, from 0.1 to 1.5mm  
10% transmittance filter  
Aperture bracket for optical bench, if necessary  
Compatible with the light sensor  
Area allowing viewing experiences.



**Laboratory Instrument India**