



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cam And Follower Mechanism

Product Code :
LBNY-0005-13000033



Description :

The unit is assembled in an anodized aluminum profile structure, with painted steel panel.

Technical Specification :

This Unit Is Mainly Composed Of:

Four cams (aluminum) of different shapes.

One flat-shaped follower (aluminum), constituting a mechanism which transforms the circular motion of the cam into a linear motion.

One eccentric (aluminum).

Two roller-shaped followers (brass) with different diameter, constituting mechanisms which transform the circular motion of the cam into the angular motion of the follower.



Laboratory Instrument India