



Email : sales@laboratoryinstrumentindia.com

Product Name :
Whitworth Quick Return Mechanism

Product Code :
LBNY-0005-13000041



Description :

This system contains bearings, flywheel and DC motor with variable speed. Easy to handle, this apparatus is preferred for its precise shape, easy to understand mechanism, application of latest technical knowledge, competitive price, low maintenance charge, sturdy construction and user friendly design. Whitworth Quick Return Mechanism is used for conversion of rotary motion into reciprocating movement at relatively slower rate than slider and crank. This apparatus is used to design equipments like mechanical actuator, screw press, shaper and power controlled saw.

Features:

- Application specific design
- Easy to handle
- Simple to understand mechanism.

Technical Specification :

Whitworth Quick Return Mechanism



Laboratory Instrument India

