



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Scotch Yoke Mechanism

**Product Code :**  
LBNY-0005-13000045



**Description :**

The main metallic elements are made of anodized aluminum.  
The unit is a bench top-unit with an anodized aluminum profile structure, and painted steel panel.  
This unit is composed of:  
A piston with a millimeter scale to measure the linear motion.  
A vertical pivot links the rotary element (crank) and the piston.  
A rotary element (crank) with a graduated disc to read the angle.

**Technical Specification :**

Scotch Yoke Mechanism



**Laboratory Instrument India**