



Email : sales@laboratoryinstrumentindia.com

Product Name :
Scotch Yoke Mechanism

Product Code :
LBNY-0005-13000045



Description :

The main metallic elements are made of anodized aluminum.
The unit is a bench top-unit with an anodized aluminum profile structure, and painted steel panel.
This unit is composed of:
A piston with a millimeter scale to measure the linear motion.
A vertical pivot links the rotary element (crank) and the piston.
A rotary element (crank) with a graduated disc to read the angle.

Technical Specification :

Scotch Yoke Mechanism



Laboratory Instrument India