



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Linear Translator

**Product Code :**  
LBNY-0004-780037



**Description :**

Linear Translator

**Technical Specification :**

Optical bench 3 accessory 06:

Linear translator:

The Linear Translator transforms a Rotary Motion Sensor into a linear motion device.

The toothed rack of the Linear Translator fits into the slot in the side of the Rotary Motion Sensor. As the Rotary Motion Sensor pulley is rotated by hand, the Rotary Motion Sensor moves along the rack. Rotary Motion Sensors not included.

Toothed nylon rack provides smooth linear motion.

Mountable on the optic bench or on rod.

For use with a rotary motion sensor.

Variable sliding distance stops.



**Laboratory Instrument India**

