




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Overdrive Unit

**Product Code :**  
LBNY-0005-13000052



#### **Description :**

Bench-top unit supported by four legs. Besides, it is provided with four brackets to be suspended in the wall, which allows the experiments with weights to be carried out more easily.

The unit is assembled in an anodized aluminum structure with a painted steel panel.

This unit is mainly composed of:

The pinions are made of aluminum and are mounted on ball bearings to reduce losses due to friction, making the practical exercises execution easy.

The shafts of the unit are made in stainless steel. The discs and gears are made in aluminum to facilitate the experiments.

2 Graduated discs located at the input and output shafts.

The epicyclic gears set, formed by:

Planet gears set.

Ring.

Sun gear.

#### **Technical Specification :**

Overdrive Unit

---

# Laboratory Instrument India

