



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Three Hinged Arch Unit

**Product Code :**  
LBNY-0005-1310008



**Description :**

This apparatus is especially suitable for measuring moving load, evenly distributed load and point load. Its left sided hinge promotes trouble free rotation and horizontal movement of its arch. Three Hinged Arch Unit is used to study unsymmetrical and symmetrical design of its structure, to conduct comparison based analysis of measured and calculated horizontal thrust during lifting load and so on. Interchangeable feature of its arch segment, application specific design and simple maintenance process are its key attributes. Side hinges located in its arch section provide required support and help to form its three hinge design. We are an eminent manufacturer and supporter of Three Hinged Arch Unit.

**Features:**

- Ergonomic structure
- Simple to handle
- Long working life
- Cost effectiveness.

**Technical Specification :**

Three Hinged Arch Unit

---

# Laboratory Instrument India

