



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Deflection Of Curved Bars Unit

**Product Code :**  
LBNY-0005-13100012



**Description :**

Bench-top unit mounted on a structure of anodized aluminum, with painted steel panel, and with regulable in heigh legs.

**Technical Specification :**

The unit basically consist of:

Structure of aluminum and steel on which the dial gauges are mounted and clamps for the test bars.

4 Test curved bars, made of steel, with different lengths and curvatures.

2 Dial gauges with range from 0-10 mm, and with 0.01 mm accuracy, to measure the deformation of the test bar subjected to a force.

Steel hook to hand the weights and it is hanging in the application point of the force.

Set of weights, composed

1 weight of 100 gr.

5 weights of 200 gr. each one.



**Laboratory Instrument India**

