




Email : sales@laboratoryinstrumentindia.com

Product Name :
Parabolic Arch Unit

Product Code :
LBNY-0005-13100013



Description :

The Unit Includes:

Parabolic arch made of tempered steel with 7 cords with hangers to hang the weights from them (uniformly distributed load or point load):

Cross section: 20 x 1.5mm.

Length: 500mm.

Height: 200mm.

Frame, of anodized aluminum, with measurement system to level the ends of the parabolic arch.

Technical Specification :

Support Systems:

Fixed end.

3 pairs of pulleys that allow to lodge the weights required for:

To cancel out the horizontal displacement and to fix the movable support.

Movable end that moves horizontally along a guide. This end can be transformed into a fixed end.

2 dial gauges, range: 0-25mm. Precision: 0.01mm.

To compensate the vertical reaction of the support.

Laboratory Instrument India

