



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hookes Law Apparatus

Product Code :
LBNY-0005-13100017



Description :

Unit mounted on an anodized aluminum frame with painted steel panel.

Technical Specification :

The Unit Includes:

2 springs with the following features:

Length: 100 mm.

Wire diameter: 1.2 mm.

External diameter: 18 mm.

1 spring with the following features:

External diameter: 18 mm.

Wire diameter: 1.1 mm.

Length: 100 mm.

These springs can be configured in series and in parallel.

Set of masses:

12 masses of 50 g.(0.5n)

2 rules to measure deformations, graduations: 1mm.

1 hook of 100g. (1n).

Laboratory Instrument India

