



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hookes Law Apparatus

**Product Code :**  
LBNY-0005-13100017



**Description :**

Unit mounted on an anodized aluminum frame with painted steel panel.

**Technical Specification :**

**The Unit Includes:**

2 springs with the following features:

Length: 100 mm.

Wire diameter: 1.2 mm.

External diameter: 18 mm.

1 spring with the following features:

External diameter: 18 mm.

Wire diameter: 1.1 mm.

Length: 100 mm.

These springs can be configured in series and in parallel.

Set of masses:

12 masses of 50 g.(0.5n)

2 rules to measure deformations, graduations: 1mm.

1 hook of 100g. (1n).

---

# Laboratory Instrument India

