



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Unit For Studying Deformation Of Trusses

**Product Code :**  
LBNY-0005-13100020



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

These bars can form several plane truss structures with different angles between bars.  
Maximum bar force: 200 N  
The unit includes: Anodized aluminum frame.

**Technical Specification :**

21 bars made of PVC, so their deformations are visible:  
1 bar of 400 mm long.  
8 bars of 300 mm long.  
1 bar of 520 mm long.  
2 bars of 150 mm long.  
4 bars of 430 mm long.  
5 bars of 259 mm long.



**Laboratory Instrument India**