



Email : sales@laboratoryinstrumentindia.com

Product Name :
Unit For Studying Deformation Of Trusses

Product Code :
LBNY-0005-13100020



www.laboratoryinstrumentindia.com

Description :

These bars can form several plane truss structures with different angles between bars.
Maximum bar force: 200 N
The unit includes: Anodized aluminum frame.

Technical Specification :

21 bars made of PVC, so their deformations are visible:
1 bar of 400 mm long.
8 bars of 300 mm long.
1 bar of 520 mm long.
2 bars of 150 mm long.
4 bars of 430 mm long.
5 bars of 259 mm long.



Laboratory Instrument India