



Email : sales@laboratoryinstrumentindia.com

Product Name :
Unit For Studying Free Vibration Of A Bar

Product Code :
LBNY-0005-13100030



Description :

Anodized aluminum frame.
Bending bar:
Material: aluminum alloy.
Section: 20 x 3 mm.
Length: 650 mm.
Fixable meter ruler for bar deflection:
Range: ± 150 mm
Electronic console, including:
Digital display for the position sensor.
Position sensor connector.
Cables and accessories, for normal operation.
Set of weights:
Material: stainless steel
Position sensor to measure the oscillation amplitude.
12 weights of 100 g. each one.

Technical Specification :

Unit For Studying Free Vibration Of A Bar

Laboratory Instrument India

