



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Brinell Hardness Testing Unit

**Product Code :**  
LBNY-0005-13100033



**Description :**

Suitable for Lab environments as well as Production from Low to Medium Test Volume  
Automatic application, holding and removal of test loads  
Selectable Dwell Times  
Best combination of quality, accuracy and cost-effectiveness  
Load Cell for Measurement of Test Forces  
Accurate and fast  
No need of handling heavy weights or cleaning the messy oils.

**Technical Specification :**

Brinell Hardness Testing Unit



**Laboratory Instrument India**