



Email : sales@laboratoryinstrumentindia.com

Product Name :
Brinell Hardness Testing Unit

Product Code :
LBNY-0005-13100033



Description :

Suitable for Lab environments as well as Production from Low to Medium Test Volume
Automatic application, holding and removal of test loads
Selectable Dwell Times
Best combination of quality, accuracy and cost-effectiveness
Load Cell for Measurement of Test Forces
Accurate and fast
No need of handling heavy weights or cleaning the messy oils.

Technical Specification :

Brinell Hardness Testing Unit



Laboratory Instrument India