




Email : sales@laboratoryinstrumentindia.com

Product Name :
LOGO! 8 PLC Trainer

Product Code :
LBNY-0002-610007



Description :

LOGO! 8 PLC Trainer

Technical Specification :

12 Inputs with switches, LEDs and 4 mm safety sockets
8 Outputs with LEDs and 4 mm safety sockets
36-pin. Centronics system connector
Power supply: 24 V DC
Suitable for A4 exercise frames or for table-top use
Dimensions: 420 x 210 x 200 mm
Weight: 2.5 kg
Including programming software

LOGO! 8 PLC Trainer, LOGO! 8 PLC Trainer Equipments, LOGO! 8 PLC Trainer Tools, LOGO! 8 PLC Trainer Tool Kits, LOGO! 8 PLC Trainer Manufacturers, LOGO! 8 PLC Trainer Suppliers from India, China, Kenya



Laboratory Instrument India

