



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wireless Pressure Sensor

Product Code :
LBNY-0004-780041



Description :

Wireless Pressure Sensor

Technical Specification :

Wireless Pressure Sensor

Bluetooth and USB connectivity with a long-lasting rechargeable battery

Logs pressure data directly onto the sensor for long-term experiments for activities that couldn't be performed using sensors previously.

Model a human lung

Includes a syringe, tubing and the connectors and locks you need for an array of pressure experiments.

Conduct weather experiments and air foil studies

Measure transpiration rates, osmosis and pressure changes inside a model cell, rates of enzyme catalyzed reactions, changes in elevation, chemical reaction rates, respiration rates etc.

No interface required! Connects directly to your devices via Bluetooth Smart.

Study Vapor Pressure vs Temperature

Pressure range : 0-400 kPA

Resolution: 0.1 kPA

Accuracy: +/- 1kPA

Connection via bluetooth

Provided with adequate adapter for connection to PC or laptop.



Laboratory Instrument India