



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computer Controlled Ventilation Trainer

**Product Code :**  
LBNY-0005-13100043



**Description :**

The designed for training in ventilation engineering. The individual components in the ventilation system are fitted with pressure measuring glands. The power consumption of the fan is indicated on the switch cabinet using a voltmeter and an ammeter. For this purpose, an inclined tube manometer with 2 measuring ranges and a digital handheld unit for the measurement of the air velocity are included, along with numerous other components. Once the trainee has built a ventilation system independently, the effects of bends, branches, regulator flaps on the airflow rate, the air velocity, as well as the pressure loss can be investigated. The components are the same as those used in ventilation engineering.

**Technical Specification :**

Computer Controlled Ventilation Trainer



**Laboratory Instrument India**