



Email : sales@laboratoryinstrumentindia.com

Product Name :
Creep Testing Unit

Product Code :
LBNY-0005-13100044



Description :

Table unit for the investigation of the creep behaviour of different materials
Tensile stress range 5...25N/mm²
Position measurement gauge 0...10mm, resolution 0.01mm
Experiments above or below room temperature possible with conditioning box
Flat specimens made of lead and PE
Specimen cross-section 2x5mm, measured length 25mm.

Technical Specification :

Specimen
Material: lead, PE
Measured length: 25mm
Tensile stress range: 5...25N/mm²
Cross-section: 2x5mm
Position measurement
Resolution 0.01mm
Range: 0...10mm.

Laboratory Instrument India

