



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifugal Casting Building Up Training Set 2

Product Code :
LBNY-0005-13100045



Description :

Aluminium cope and drag box, screwable
Pouring ladle with 5 kg Woods metal (melting point 70C)
Worksheets for instructor and students
Pattern plate with pattern of epoxy-resin in one piece
Core die, core iron
15 kg moulding sand in a separate container
Moulding tools
Thermometer.

Technical Specification :

Centrifugal Casting Building Up Training Set 2



Laboratory Instrument India