



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Centrifugal Casting Building Up Training Set 2

**Product Code :**  
LBNY-0005-13100045



**Description :**

Aluminium cope and drag box, screwable  
Pouring ladle with 5 kg Woods metal (melting point 70C)  
Worksheets for instructor and students  
Pattern plate with pattern of epoxy-resin in one piece  
Core die, core iron  
15 kg moulding sand in a separate container  
Moulding tools  
Thermometer.

**Technical Specification :**

Centrifugal Casting Building Up Training Set 2



**Laboratory Instrument India**