



Email : sales@laboratoryinstrumentindia.com

Product Name :
Deadweight Tester

Product Code :
LBNY-0005-1320003



Description :

A simple deadweight tester with weights and a bourdon tube pressure gauge with a clear front allowing the mechanism to be viewed. The unit has an optional pressure transducer available allowing the concept of electronic representation of pressure to be investigated. Unique design allows standard Newton weights to be used; hence reinforcing the student understanding of pressure being Force / Unit Area.

Technical Specification :

Deadweight Tester



Laboratory Instrument India