



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Compressible Flow Pipe Friction Module

**Product Code :**  
LBNY-0005-1320004



**Description :**

The pressures, temperatures and air flow rate are recorded by a combination of instrumentation on the Compressible Flow Range base unit and the optional module. A digital hand held manometer is supplied. A series of four straight tubes of different diameters with end pressure tapings allow pressure losses in a straight pipe to be investigated at a range of Reynolds numbers. Bends, sudden enlargements and contractions are included to allow investigation of pressure loss and recovery. The unit is driven by an ejector (jet pump) allowing investigation of entrainment ratio and ejector performance.

**Technical Specification :**

Compressible Flow Pipe Friction Module



**Laboratory Instrument India**