




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Reaction Turbine Module

**Product Code :**  
LBNY-0005-1320005



**Description :**

Inlet and outlet air pressures, temperatures and air flow rate, turbine torque and speed are recorded by a combination of instrumentation on the Compressible Flow Range base unit and the optional module.  
Experimental Reaction Turbine Module:- Inlet and outlet air temperatures are recorded allowing the temperature drop due to work output to be measured. A single stage, radial flow , two jet reaction turbine with a throttle valve and belt brake dynamometer.

**Technical Specification :**

Reaction Turbine Module



**Laboratory Instrument India**