



Email : sales@laboratoryinstrumentindia.com

Product Name :
Nozzle Performance Test Module

Product Code :
LBNY-0005-13200010



Description :

A pressure regulator, throttle valve and back pressure valve allow the air flow rate, inlet and discharge (or back) pressure to be varied. A series of convergent and convergent-divergent nozzles may be installed in one of two locations in a high pressure measuring chamber. Inlet and outlet air pressures, temperatures and air flow rate are recorded by a combination of instrumentation on the Compressible Flow Range base unit and the optional module. Standard unit includes convergent & divergent ducts designed to produce Mach 1.0 at the throat and supersonic velocities downstream.

Technical Specification :

Nozzle Performance Test Module



Laboratory Instrument India