



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Impulse Turbine Module

**Product Code :**  
LBNY-0005-13200011



**Description :**

Inlet and outlet air pressures, temperatures and air flow rate, turbine torque and speed are recorded by a combination of instrumentation on the Compressible Flow Range base unit and the optional module. Inlet and outlet air temperatures are recorded allowing the temperature drop due to work output to be measured. An impulse turbine with 4 separate nozzles and control valves, a throttle valve and belt brake dynamometer.

**Technical Specification :**

Impulse Turbine Module



**Laboratory Instrument India**