



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Rubber In Shear Apparatus

**Product Code :**  
LBNY-0005-1330008



**Description :**

The unit has dial gauge attached to the upper end of the rubber block to measure the amount of deformation. The unit consists of a frame that is attached to wall. To produce shear strain in the rubber block, load is applied to metallic strip via a hanger. A rubber block is affixed to its surface. The rubber block has a metallic strip fastened to its surface.

This apparatus is designed to study the shear strain in a rubber block due to the shear stress applied on it. When a force acts parallel to the surface of any object, it produces shear strain in it.

**Technical Specification :**

Rubber In Shear Apparatus



**Laboratory Instrument India**